

FACE MASK

DISPOSABLE

Protective | Comfort Fit | Safe and Secure



*FOR REFERENCE ONLY



Professional Filter Material



3-Layer Protection



High Efficiency & Low Resistance



Comfortable & Breathable

3-LAYER OF PROTECTION

BFE \geq 95%



- Pleated, soft, and skin-friendly.
- Flexible nose wire and elastic earloop.
- Effective against bacteria, aerosols, pollen, and dust.
- Waterproof and odor-free, comfortable, and breathable.
- Applying for FDA approval and CE mark.

PRODUCT DETAILS

Product Name	Disposable Face Mask
Composition	Non-woven fabric, meltblown fabric, ear loop, and nose wire
Specifications	Bacterial Filtration Efficiency \geq 95% Breathing Resistance \leq 49Pa/Cm ² (at 8L / Min)
Usage	Used for daily respiratory protection and blocking of pollutants
Shelf Life	1 year after production date
Manufacturer	Shenzhen Hytera Communications Corporation Limited
Address	4th floor, Plant 2, Hytera Hi-tech Park No.3, Baolong 4th Rd., Baolong Industrial Area, Longgang District, Shenzhen, Guangdong, P.R. China
Tel	+86 755 89788999

How to put on the mask



- 1 Place the mask over your face with the colored side facing outwards.
- 2 Place the ear loops around your ears.
- 3 Pinch the nose wire over your nose bridge.
- 4 Pull the bottom of the mask over your mouth and chin.

CHILD'S Disposable Face Mask

Protective | Comfort Fit | Safe and Secure



* FOR REFERENCE ONLY



Professional Filter Material



3-Layer Protection



High Efficiency & Low Resistance



Comfortable & Breathable

PRODUCT DETAILS

Product Name	Child's Disposable Face Mask
Composition	Non-woven fabric, meltblown fabric, ear loop, and nose wire
Specifications	Bacterial Filtration Efficiency $\geq 95\%$ Breathing Resistance $\leq 49\text{Pa}/\text{Cm}^2$ (at 8L / Min)
Usage	Used for daily respiratory protection and blocking of pollutants
Shelf Life	1 year after production date
Manufacturer	Shenzhen Hytera Communications Corporation Limited
Address	4th floor, Plant 2, Hytera Hi-tech Park No.3, Baolong 4th Rd., Baolong Industrial Area, Longgang District, Shenzhen, Guangdong, P.R. China
Tel	+86 755 89788999

How to put on the mask

- 1 Place the mask over your face with the colored side facing outwards.
- 2 Place the ear loops around your ears.
- 3 Pinch the nose wire over your nose bridge.
- 4 Pull the bottom of the mask over your mouth and chin.



3-LAYER OF PROTECTION

BFE $\geq 95\%$



- Pleated, soft, and skin-friendly.
- Flexible nose wire and elastic earloop.
- Effective against bacteria, aerosols, pollen, and dust.
- Waterproof and odor-free, comfortable, and breathable.
- Applying for FDA approval and CE mark.





默认发货为无单独塑料薄膜袋包装，如需单独包装需备注

Quantity(PCS)	Unit price (USD/PCS)	(USD/PCS)	Price term
10K-100M (Including 10K)	0.60		FOB Shenzhen
100M – 500M (Including 100M)	0.50		FOB Shenzhen
More than 500M (Including 500M)	0.45		FOB Shenzhen